ALIAXIN® SV


Product Specifications

1-844-226-8277
info@medsupplysolutions.com
medsupplysolutions.com

Aliaxin[®]SV

SUPERIOR VOLUME


Cross-Linked
Hyaluronic Acid


hydroACTION ●○○○○○
liftACTION ●●●●●●
CROSS-LINKING ■■■■■■


1


2


3

4

5

IT

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COMPOSITION

Hyaluronic acid * 25 mg / mL

*In the form of cross-linked sodium hyaluronate

Physiological pH phosphate buffer pH QS

DESCRIPTION

Aliaxin® SV Superior Volume is a viscoelastic hyaluronic acid gel, sterile, apyrogenic and of non-animal origin.

INDICATIONS

Aliaxin® SV Superior Volume is indicated for:

- volume restoration and the filling of cutaneous depressions and wrinkles
- facial reconstruction of volumes in the event of light-to-moderate lipodystrophies, i.e. due to HIV treatment

Aliaxin® SV Superior Volume is indicated for the augmentation of cheeks, cheekbones and facial contours, as well as the filling of very deep wrinkles.

STERILISATION

The content of the Aliaxin® SV Superior Volume syringes is sterilised using moist heat. The needles are sterilised using radiation

INTENDED USE

Aliaxin® SV Superior Volume is a filling gel that is injected into the subcutaneous layer in order to restore volume in the required proportions.

CONTRAINDICATIONS

Aliaxin® SV Superior Volume must not be used:

- for injections other than subcutaneous ones, and more specifically it absolutely must not be injected into blood vessels. Do not inject intravascularly.
- in the periorbital and glabellar area.

Do not inject in cases of:

- Known hypersensitivity to any of the components,
- A history of severe allergies or anaphylactic shock.
- In combination with skin peeling treatment, laser treatment, ultrasound-based treatment or a technique based on an active dermal response.
- Inflammatory and/or infectious cutaneous problems on or near the area to be treated (acne, herpes, etc.) or tumours

Do not inject into patients:

- Suffering from autoimmune diseases,
- Suffering from epilepsy which is not controlled by treatment,
- That tend to develop hypertrophic scars
- Suffering from porphyria
- That are pregnant or breast-feeding
- That are under 18 years of age.

PRECAUTIONS FOR USE

Before starting treatment, patients should be questioned about their medical history, and be informed of the treatment's predictable results and its potential adverse effects.

It is not advisable to inject Aliaxin® SV Superior Volume products into patients that have taken vitamin E or vitamin C in large doses, aspirin, anti-inflammatories or anti-coagulants.

This device is intended for single use only. Sterility is not guaranteed in the event of re-use.

Do not use the product if the box is damaged or if the expiration date has passed.

It is not advisable to inject into an area which has already been treated with a semi-permanent or permanent implant.

POST-TREATMENT CARE

After treatment, patients must be advised:

- Not to apply make-up in the 12 hours following the injection
- Not to expose themselves in the two weeks following the treatment to extreme temperatures (intense cold, sauna, hammam, etc.) or sunlight, nor to take sunbed sessions in the week following the treatment.

Do not apply pressure to the treated area.

DOSAGE

Aliaxin® SV Superior Volume may only be injected by practitioners who are trained in the techniques of injecting filling products and who are authorised to perform injections.

Practitioners must inject the amount depending on the desired effect, but without over-correcting: they must carry out at the most 100% of the desired correction, without going further.

Aliaxin®SV Superior Volume is a biodegradable gel that is absorbed over time and may require several injections before a satisfactory result is obtained.

The practitioner may suggest subsequent sessions to maintain the desired effect. The frequency of these sessions is determined by the results of previous injections.

The average duration depends on several factors: the type of patient's skin, the severity of the wrinkle to be corrected, type of injection (superficial, medium or deep dermis) and the injected volume.

MODE OF ADMINISTRATION

Before injection:

- Disinfect the treatment area.
- Remove the air by gently pushing down on the plunger to make a small amount of the product appear at the tip of the needle.
- Do not mix with other products before injection.

Do not bring together with quaternary ammonium salts due to a risk of precipitate. In particular, do not bring together with equipment sterilised with benzalkonium chloride.

Injection:

- It is recommended to perform the injection using the 27G needles provided in the box.
- Aliaxin®SV Superior Volume is injected evenly on the areas to be treated into the subcutaneous layer
- An injection which is too superficial or one whose injected volume is excessive may cause bumps on the treated area.

After injection:

- After the injection, it is recommended to massage the treated area to enable the product to spread evenly.
- During the injection, if a sensation of obstruction or pressure is felt, the injection must be stopped and the needle must be changed.
- Always dispose of the syringe, needle and unused product immediately after the treatment and do not use them again. This is to prevent any contamination of the remaining product and to avoid the associated risks, especially risks of infection.

These items must be disposed of in accordance with applicable medical practice as well as with institutional or local and national rules.

SIDE EFFECTS

The patient must be informed that:

- symptoms of inflammation on the implant site appear soon after the injection or with a delayed onset.
- in general, the injection of a medical device carries the risk of infection.

Reactions that may occur due to the injection are as follows:

- Oedema
- Haematoma,
- Erythema,
- Itchiness,
- Swelling,
- Pain or tenderness at the implant site,
- Induration or nodules at the injection point, perhaps leading to granuloma.
- Cutaneous necrosis
- Immediate hypersensitivity, perhaps leading to anaphylactic shock
- Dyschromia at the injection site

If inflammatory reactions persist beyond one week or if any other side effect occurs, the patient should inform the practitioner as soon as possible. The practitioner should then carry out the appropriate treatment.

Adverse effects should be reported to manufacturer and/or local distributor of the product.

STORAGE CONDITION

The expiration date is indicated on the packaging.

Store between 2°C and 25°C, away from light. Do not freeze.

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DESCRIPTION

Aliaxin® EV Essential Volume is a reabsorbable medical device (physiological, apyrogenous sterile gel) to be used as a filler for the correction of deep facial cutaneous sagging and to increase volume. The principle component is cross-linked Hyaluronic Acid of non-animal origin, produced by bacterial fermentation.

COMPOSITION

Cross-linked hyaluronic acid.....25 mg/g
Phosphate buffered saline, Water for injectable solutions.....1 g

PACKAGE

Aliaxin® EV Essential Volume is available in a 2ml package. Each package contains:

- package insert leaflet
- 2 sealed blisters, each containing 1 sterile pre-filled mono-use syringe of 1 ml
- Adhesive labels to stick to the patient's medical record for product traceability
- 4 sterile needles (27G X 13mm TSK C€0123 - 27G X 19mm TSK C€0123)

MODALITY OF ACTION

Aliaxin® EV Essential Volume must be injected into the sub-cutaneous tissue for the supplementation of the extra-cellular root and the volumetric increase of the tissue, that is, for the treatment of cutaneous depressions caused by wrinkles, scarring and volume loss.

AREAS OF USE

Aliaxin® EV Essential Volume is an injectable medical device conforming to the Dir. 93/42/CEE MOD, indicated for the treatment of imperfections in the following areas of the face:

- nasolabial folds
 - post-traumatic and acne scars
 - areas of the face which require enrichment of the facial tissue (cheeks, chin, cheekbones) with temporary volumetric increase through growth of the soft tissues
- The outcome of the treatment depends on the type of skin and the nature of the imperfection: the less evident the imperfection, the better the result.

INSTRUCTIONS FOR USE

Before any form of intervention with Aliaxin® EV Essential Volume, the doctor must carry out an adequate case history and a comprehensive assessment of the conditions of the patient, in order to guarantee the total absence of contraindications to the implant.

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Local anaesthetic may be used in order to guarantee the necessary comfort to the patient.

The doctor must inform the patient in advance of the modalities of the intervention, its nature, warnings, precautions and possible individual outcomes, potential adverse responses, the expected duration of the treatment and the eventuality of a supplementary intervention for the maintenance and/or definition in detail of the achieved result.

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Extract the syringe from the blister pack, unscrew the top and assemble the needle together with the protection. The latter is only to be removed right before the intervention. Aliaxin® EV Essential Volume is administered with a sterile needle (27G X 13mm - 27G X 19mm) in compliance with the Luer-Lock standards with normalised connections, which is included in the pack.

The injection must be carried out in correspondence to the sub-cutaneous superficial layer; the procedure is, however, at the discretion of the doctor and depends on the correction desired and the method of implantation adopted.

At the end of the procedure it is advisable to carry out a delicate massage of the treated area in order to allow for better distribution of the product.

DOSAGE AND ADMINISTRATION

The volumetric grading printed on the syringe has an indicative value: it is left to the doctor's discretion to choose the dosage to be used for each individual case. The intervals at which the treatment should be repeated depend on various factors, regarding both the physiology of the patient (type of skin, individual metabolism, anatomy, age) and their lifestyle. Another element to consider is relative to the implant techniques used. For maintenance of the obtained results it is advisable to repeat the treatment periodically, every 8-12 months.

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Aliaxin® EV Essential Volume must not be used on patients affected by:

- infective or inflammatory processes near the area of intervention
- ongoing cutaneous allergies
- known hypersensitivity to cheloids
- allergy to the components
- immune system illnesses
- serious pathological skin conditions
- problems with coagulation or anti-coagulating therapies in progress

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Aliaxin® EV Essential Volume may not be implemented in areas where there are mammary, tendon, bone or muscular implants present.

Following the treatment, and until the complete absorption of swelling and reddening, the implanted areas must not be exposed to make-up, excessive heat (sun, UV tanning sessions or laser) or intense cold.

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Aliaxin® EV Essential Volume is an intra cutaneous injection and must not be injected into blood vessels.

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Do not use in pregnancy.

Aliaxin® EV Essential Volume is available in single dose sterile packages.

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Do not mix with other injectable substances and do not carry out other implants together with Aliaxin® EV Essential Volume.

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Fill in the adhesive label contained in the box and stick it to the patient's medical records kept at the doctor's office.

Keep out of the reach of children.

STORAGE

Aliaxin® EV Essential Volume must be stored at a temperature between 2° and 28°C. Do not freeze. Do not expose to heat.

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DESCRIPTION

Aliaxin® LV Lips Volume is a reabsorbable medical device (physiological, apyrogenous sterile gel) to be used as a filler for the correction of medium and deep facial cutaneous sagging and to increase the volume and contours of the lips. The principle component is cross-linked Hyaluronic Acid of non-animal origin, produced by bacterial fermentation.

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PACKAGE

Aliaxin® LV Lips Volume is available in a 2ml package. Each package contains:

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MODALITY OF ACTION

Aliaxin® LV Lips Volume must be injected into the cutaneous tissue for the supplementation of the extra-cellular root and the intra-cutaneous volumetric increase, that is, for the treatment of cutaneous depressions caused by wrinkles and scarring.

AREAS OF USE

Aliaxin® LV Lips Volume is an injectable medical device designed to restore physiological facial volume and indicated for the treatment of imperfections in the following areas of the face:

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INSTRUCTIONS FOR USE

Before any form of intervention with Aliaxin® LV Lips Volume, the doctor must carry out an adequate case history and a comprehensive assessment of the conditions of the patient, in order to guarantee the total absence of contraindications to the implant.

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Extract the syringe from the blister pack, unscrew the top and assemble the needle together with the protection. The latter is only to be removed right before the intervention. Aliaxin® LV Lips Volume is administered with a sterile needle (27G X 13mm - 27G X 19mm) in compliance with the Luer-Lock standards with normalised connections, which is included in the pack.

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Aliaxin® LV Lips Volume must not be used on patients affected by:

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Aliaxin® LV Lips Volume may not be implemented in areas where there are mammary, tendon, bone or muscular implants present.

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Aliaxin® LV Lips Volume is an intra cutaneous injection and must not be injected into blood vessels.

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Do not use in pregnancy.

Aliaxin® LV Lips Volume is available in single dose sterile packages.

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Do not mix with other injectable substances and do not carry out other implants together with Aliaxin® LV Lips Volume.

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Keep out of the reach of children.

STORAGE

Aliaxin® LV Lips Volume must be stored at a temperature between 2° and 28°C. Do not freeze. Do not expose to heat.

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DESCRIPTION

Aliaxin® FL Lips is a reabsorbable medical device (physiological, apyrogenous sterile gel) to be used as a filler for the correction of medium and deep facial cutaneous sagging and to increase the volume and contours of the lips. The principle component is cross-linked Hyaluronic Acid of non-animal origin, produced by bacterial fermentation.

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PACKAGE

Aliaxin® FL Lips is available in a 2ml package. Each package contains:

- package insert leaflet
- 2 sealed blisters, each containing 1 sterile pre-filled mono-use syringe of 1 ml
- adhesive labels to stick to the patient's medical record for product traceability
- 4 sterile needles (30G X 13mm TSK C € 0123)

MODALITY OF ACTION

Aliaxin® FL Lips must be injected into the cutaneous tissue for the supplementation of the extra-cellular root and the intra-cutaneous volumetric increase, that is, for the treatment of cutaneous depressions caused by wrinkles and scarring.

AREAS OF USE

Aliaxin® FL Lips is an injectable medical device designed to restore physiological facial volume and indicated for the treatment of imperfections in the following areas of the face:

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INSTRUCTIONS FOR USE

Before any form of intervention with Aliaxin® FL Lips, the doctor must carry out an adequate case history and a comprehensive assessment of the conditions of the patient, in order to guarantee the total absence of contraindications to the implant. The areas to be treated must be identified and evaluated bearing in mind the criteria relative to lines and symmetries to follow.

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Extract the syringe from the blister pack, unscrew the top and assemble the needle together with the protection. The latter is only to be removed right before the intervention.

Aliaxin® FL Lips is administered with a sterile needle (30G X 13mm) in compliance with the Luer-Lock standards with normalised connections, which is included in the pack.

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Aliaxin® FL Lips must not be used on patients affected by:

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Aliaxin® FL Lips may not be implemented in areas where there are mammary, tendon, bone or muscular implants present.

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Aliaxin® FL Lips is an intra cutaneous injection and must not be injected into blood vessels.

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Do not use in pregnancy.

Aliaxin® FL Lips is available in single dose sterile packages.

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Do not mix with other injectable substances and do not carry out other implants together with Aliaxin® FL Lips.

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Keep out of the reach of children.

STORAGE

Aliaxin® FL Lips must be stored at a temperature between 2° and 28°C.

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The product is only to be used by medical staff as described in this information leaflet. Do not use for applications different than those indicated in this leaflet.

DESCRIPTION

Aliaxin® GP Global Performance is a reabsorbable medical device (physiological, apyrogenous sterile gel) to be used as a filler for the correction of medium and deep facial cutaneous sagging and to increase the volume and contours of the lips. The principle component is cross-linked Hyaluronic Acid of non-animal origin, produced by bacterial fermentation.

COMPOSITION

Cross-linked hyaluronic acid.....25 mg/g
Phosphate buffer, Water for injectable solutions 1 g

PACKAGE

Aliaxin® GP Global Performance is available in a 2ml package. Each package contains:

- package insert leaflet
- 2 sealed blisters, each containing 1 sterile pre-filled mono-use syringe of 1 ml
- Adhesive labels to stick to the patient's medical record for product traceability
- 4 sterile needles (27G X 13mm TSK C € 0123 - 27G X 19mm TSK C € 0123)

MODALITY OF ACTION

Aliaxin® GP Global Performance must be injected into the cutaneous tissue for the supplementation of the extra-cellular root and the intra-cutaneous volumetric increase, that is, for the treatment of cutaneous depressions caused by wrinkles and scarring.

AREAS OF USE

Aliaxin® GP Global Performance is an injectable medical device designed to restore physiological facial volume and indicated for the treatment of imperfections in the following areas of the face:

- nasolabial folds
- wrinkles around the mouth
- facial and acne scars
- glabellar wrinkles
- wrinkles on the forehead
- areas of the face which require enrichment of the facial tissue (cheeks, chin, cheekbones, lips) with temporary volumetric increase through growth of the soft tissues

The outcome of the treatment depends on the type of skin and the nature of the imperfection: the less evident the imperfection, the better the result.

INSTRUCTIONS FOR USE

Before any form of intervention with Aliaxin® GP Global Performance, the doctor must carry out an adequate case history and a comprehensive assessment of the conditions of the patient, in order to guarantee the total absence of contraindications to the implant.

The areas to be treated must be identified and evaluated bearing in mind the criteria relative to lines and symmetries to follow.

Local anaesthetic may be used for treatment to the lips, in order to guarantee the necessary comfort to the patient.

The doctor must inform the patient in advance of the modalities of the intervention, its nature, warnings, precautions and possible individual outcomes, potential adverse responses, the expected duration of the implant and the eventuality of a supplementary intervention for the maintenance and/or definition in detail of the achieved result.

The area of intervention must be cleansed with antiseptic solution before proceeding to the implant.

Extract the syringe from the blister pack, unscrew the top and assemble the needle together with the protection. The latter is only to be removed right before the intervention. Aliaxin® GP Global Performance is administered with a sterile needle (27G X 13mm - 27G X 19mm) in compliance with the Luer-Lock standards with normalised connections, which is included in the pack.

The implantation must be carried out in correspondence to the medium/deep derma; the procedure is, however, at the discretion of the doctor and depends on the correction desired and the method of implantation adopted.

At the end of the procedure it is advisable to carry out a delicate massage of the treated area in order to allow for better distribution of the implant.

DOSAGE AND ADMINISTRATION

The volumetric grading printed on the syringe has an indicative value: it is left to the doctor's discretion to choose the dosage to be used for each individual case. The intervals at which the treatment should be repeated depend on various factors, regarding both the physiology of the patient (type of skin, individual metabolism, anatomy, age) and their lifestyle. Another element to consider is relative to the implant techniques used. For maintenance of the obtained results it is advisable to repeat the treatment periodically, every 6-12 months.

PRECAUTIONS

During the treatment the usual precautions in the case of percutaneous procedures must be adopted. The risks are the common ones of infection relative to the type of treatment.

Aliaxin® GP Global Performance must not be used on patients affected by:

- infective or inflammatory processes near the area of intervention
- known hypersensitivity to cheloids
- allergy to the components
- immune system illnesses
- serious pathological skin conditions
- problems with coagulation or anti-coagulating therapies in progress

Close to treatment the patient should avoid consumption of substances which act on blood fluidity (aspirin, NSAID, Vitamin E), in order to minimise the possibility of bruising or bleeding of the injected areas.

Aliaxin® GP Global Performance may not be implemented in areas where there are mammary, tendon, bone or muscular implants present.

Following the treatment, and until the complete absorption of swelling and reddening, the implanted areas must not be exposed to excessive heat (sun, UV tanning sessions or laser) or intense cold.

After use, needles and syringes must be disposed of according to the procedures indicated for hospital waste.

ADVERSE REACTIONS, WARNINGS AND CONTRAINDICATIONS

As with any percutaneous injection, there can be incidents of:

- inflammatory reactions (reddening, edema, etc), sometimes associated with itching and pain to the touch
- hematoma
- hardening or nodules at the point of injection
- colouring or discolouring of the skin at the point of injection

These reactions usually disappear within a few days. If they persist for more than a week the patient should contact their doctor immediately.

Aliaxin® GP Global Performance is an intra cutaneous injection and must not be injected into blood vessels.

There are no reports of problems due to excessive doses or interaction with other medicines.

Do not use in pregnancy.

Aliaxin® GP Global Performance is available in single dose sterile packages.

The product must not be reused for future applications or on other patients.

Once open the product must be used immediately.

Any remaining product must be disposed of.

Do not use the product if the packaging is damaged.

Do not mix with other injectable substances and do not carry out other implants together with Aliaxin® GP Global Performance.

The implant must be carried out in an aseptic environment observing the appropriate techniques.

Fill in the adhesive label contained in the box and stick it to the patient's medical records kept at the doctor's surgery.

Keep out of the reach of children.

STORAGE

Aliaxin® GP Global Performance must be stored at a temperature between 2° and 28°C. Do not freeze. Do not expose to heat.

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
SHAPE & RESTORE


Cross-Linked
Hyaluronic Acid


hydroACTION ●●●●○
liftACTION ●●●○○
CROSS-LINKING ■■■□□

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


EN

ASSEMBLING THE NEEDLE ON THE SYRINGE

For optimal handling of Aliaxin® SR Shape & Restore products, it is important to assemble the needle on the syringe in accordance with the 5 steps shown in the diagram (see figures 1 to 5).

1

Firmly hold the glass syringe and the luer-lock between the thumb and the index fingers. Grasp the protective cap and unscrew it.

2

Remove needle cap, then firmly insert the needle thread on the end of the syringe.

3

Gently screw the needle in a clockwise direction. Continue to screw until the thread comes into contact with the body of the syringe.

4

Check that the needle is correctly fitted onto the thread. Failure to follow these precautions could lead to risks of the needle dropping out and/or leaks from the luer-lock.

5

With the syringe held in one hand and the cap in the other, remove the needle protective device by pulling firmly.

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DESCRIPTION

Aliaxin® SR Shape & Restore is a reabsorbable medical device (physiological, apyrogenous sterile gel) to be used as a filler for the correction of medium and deep facial cutaneous sagging and to increase the volume and contours of the lips. The principle component is cross-linked Hyaluronic Acid of non-animal origin, produced by bacterial fermentation.

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PACKAGE

Aliaxin® SR Shape & Restore is available in a 2ml package. Each package contains:

- package insert leaflet
- 2 sealed blisters, each containing 1 sterile pre-filled mono-use syringe of 1 ml
- adhesive labels to stick to the patient's medical record for product traceability
- 4 sterile needles (30G X 13mm TSK C €0123 - 27G X 19mm TSK C €0123)

MODALITY OF ACTION

Aliaxin® SR Shape & Restore must be injected into the cutaneous tissue for the supplementation of the extra-cellular root and the intra-cutaneous volumetric increase, that is, for the treatment of cutaneous depressions caused by wrinkles and scarring.

AREAS OF USE

Aliaxin® SR Shape & Restore is an injectable medical device designed to restore physiological facial volume and indicated for the treatment of imperfections in the following areas of the face:

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INSTRUCTIONS FOR USE

Before any form of intervention with Aliaxin® SR Shape & Restore, the doctor must carry out an adequate case history and a comprehensive assessment of the conditions of the patient, in order to guarantee the total absence of contraindications to the implant.

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Extract the syringe from the blister pack, unscrew the top and assemble the needle together with the protection. The latter is only to be removed right before the intervention. Aliaxin® SR Shape & Restore is administered with a sterile needle (30G X 13mm - 27G X 19mm) in compliance with the Luer-Lock standards with normalised connections, which is included in the pack.

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Aliaxin® SR Shape & Restore must not be used on patients affected by:

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Aliaxin® SR Shape & Restore may not be implemented in areas where there are mammary, tendon, bone or muscular implants present.

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Aliaxin® SR Shape & Restore is an intra cutaneous injection and must not be injected into blood vessels.

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Do not use in pregnancy.

Aliaxin® SR Shape & Restore is available in single dose sterile packages.

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Do not mix with other injectable substances and do not carry out other implants together with Aliaxin® SR Shape & Restore.

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Keep out of the reach of children.

STORAGE

Aliaxin® SR Shape & Restore must be stored at a temperature between 2° and 28°C.

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 Rose Pharma S.A.
Via San Gottardo, 10
6900 Lugano - Switzerland

Rev. 5 del 03/2019

This document is a photographic image of the product insert provided for your convenience by Med Supply Solutions. Med Supply Solutions does not claim ownership of the content contained within this product insert. The intellectual property rights, including but not limited to copyrights, trademarks, and any other applicable rights, pertaining to the content of this product insert, remain the sole property of the original manufacturer. Med Supply Solutions accepts no liability for any inaccuracies, omissions, or discrepancies that may be present in this document. Your use of this product insert is subject to our terms and conditions, which can be found on our website. By accessing and utilizing this document, you acknowledge and accept the terms and conditions outlined therein.

1-844-226-8277
info@medsupplysolutions.com
medsupplysolutions.com

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IBSA Farmaceutici Italia Srl