STYLAGE® XXL
BI-SOFT®
Product Specifications

1-844-226-8277
info@medsupplysolutions.com
medsupplysolutions.com

STYLAGE®
(P) - Like TECHNOLOGY

XL, XXL
XL Lidocaïne

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DESCRIPTION

The products are cross-linked hyaluronic acid gels of non-animal origin, incorporating an antioxidant (mannitol). They are sterile and non-pyrogenic, with physiological pH and osmolality.

The gels are packaged in a pre-filled syringe and sterilised using autoclave. These products are for single use. Each box contains one instruction leaflet and labels with the batch number – one of which must be given to the patient while the other is attached to the doctor's patient file. Each box also contains:

Stylage® XL	2 x 1 ml pre-filled syringes, 4 x single-use sterile needles
Stylage® XL Lidocaïne	2 x 1 ml pre-filled syringes, 4 x single-use sterile needles
Stylage® XXL	2x 1 ml pre-filled syringes, 3 x single-use sterile needles, 1 x single-use sterile cannula or 1 x 2.2 ml pre-filled syringe, 3 x single-use sterile needles, 1 x single-use sterile cannula

INDICATIONS

Stylage® XL Stylage® XXL	A hyaluronic acid injectable gel indicated for the restoration or augmentation of facial volume by injection into the deep dermis or subcutaneously.
Stylage® XL Lidocaïne	A hyaluronic acid injectable gel indicated for the restoration or augmentation of facial volume by injection into the deep dermis or subcutaneously. Lidocaine hydrochloride is intended to reduce the pain associated with the injection.

MODE OF ACTION

The products must be injected into the area to be treated, where they will generate volume. The products will then be resorbed slowly over time. The duration of the products depends on the patient's skin type and the depth of injection. Thus, based on the characteristics of the area treated, the patient, and the depth of the injection, one to two injection sessions are necessary for optimal treatment of the aesthetic defect. Through regular touch-up sessions, it will be possible to make the desired correction more durable.

CONTRAINDICATIONS

The products must not be used:

- In combination with a peeling, laser treatment, or dermabrasion. The practitioner shall decide upon the necessary waiting period after these treatments before the injection;
- In areas where the skin is thin (e.g. the forehead and the periorbital region, including the eyelids, the under-eye shadows, the crow's feet), in areas where the vascular system is exposed (e.g. the glabellar region) and in the lips;
- In patients with a tendency to develop hypertrophic scars;
- In patients with a known hypersensitivity to one of the ingredients;
- In women who are pregnant or breastfeeding, or in children;
- In areas, or near areas, with inflammatory or infectious skin problems (acne, herpes, etc.).

Furthermore, **Stylage® XL Lidocaïne** must not be used:

- In patients with a known hypersensitivity to amide-type local anaesthetics;
- In patients with porphyria.

PRECAUTIONS FOR USE

- It is strongly recommended to obtain the patient's informed consent prior to treatment.
- If the patient has a history of herpes, there is a risk that needle punctures may trigger a new episode of herpes.
- For patients with a history of, or who are currently suffering from an autoimmune disease, the practitioner must decide on the indication on a case-by-case basis, based on the nature of the illness and the associated treatments. The practitioner must also ensure special monitoring of these patients, which includes offering a double test before injection, and he/she must not perform the injection if the illness is progressive.

- Patients with a history of streptococcal disease, such as recurrent sore throat or acute rheumatic fever, must undergo a double test before any injection. The injection is not recommended in cases of acute rheumatic fever with cardiac localisation.
 - Combining **Stylage® XL Lidocaïne** with certain medications that reduce or inhibit hepatic metabolism is not recommended.
 - In case of haemostatic disorder or anticoagulant treatment, the risk of haematoma is increased.
 - The intake of aspirin, NSAIDs, platelet aggregation inhibitors, anticoagulants, or vitamin C during the week preceding the injection should be avoided.
 - The products must not be injected into the blood vessels or muscles. Any accidental intravascular injection may cause vascular occlusion, which can lead to rare but serious complications such as vision disorders, blindness, necrosis of the skin and/or underlying tissues, depending on the area injected.
 - The products must not be injected into the nerves. Any accidental nerve damage may cause transient paresthesia.
 - It is not advised to inject implants from the **Stylage®** product line into an area that has already been corrected using implants from a product line other than **Stylage®** as no clinical data is available.
 - It is advised not to mix it with another product.
 - Only the needles/cannula provided with the products should be used for the injection, as the combination of these two devices has been validated.
 - For the injection of the **Stylage® XL** and **Stylage® XL Lidocaïne** products, the practitioner will choose the needle of a suitable length for the area to be treated, the type of correction desired and the injection technique chosen.
 - The use of cannulas with foam tips (sterile, CE marked, with a conical connection hub) is possible. In this case, the practitioner will need to:
 - . select the appropriate cannula size (gauge and length) for the correction to be made;
 - . comply strictly with aseptic conditions (particularly pre-hole).
 - Inject slowly to avoid overcorrections.
 - In case of increased pain during injection, stop the procedure and withdraw the needle/cannula.
 - Do not use if the packaging has been damaged (syringe, blister pack, lid).
 - Inject immediately after opening.
 - Discard the syringe and the residual product when the injection is complete, and dispose of the needles/the cannula in an appropriate container, in accordance with the current national regulations.
 - The **Stylage®** implants are for single use. Do not reuse.
- An implant must be used only for one single patient to prevent any risk of cross-contamination.
- After opening, the product can never be re-sterilised, even if the implant has not been injected.
 - The product **Stylage® XL Lidocaïne** contains an active ingredient, lidocaine hydrochloride, which can induce a positive result in tests performed in anti-doping controls.

DOSAGE - METHOD OF ADMINISTRATION

- The treatment should be carried out by a legally qualified practitioner, trained in implant injection techniques and with good knowledge of the anatomy of the face and the injection planes. The practitioner should also consider the presence of lidocaine hydrochloride in the product **Stylage® XL Lidocaïne**. The practitioner chooses the product to inject based on the anatomy of the area and the desired effect.
- Before carrying out the treatment, it is essential:
 - . to thoroughly confirm the patient's medical history with her/him;
 - . to explain the indications and expected results of the products to the patient;
 - . to explain the contraindications, precautions for use, and potential adverse effects related to the treatment with the patient as well as the "advice to patients".
- Before injection:
 - . disinfect the area to be treated with an appropriate antiseptic solution;
 - . verify the integrity of the needles;
 - . verify that the gel is not cloudy.
- In the event of refrigerated storage, bring the product to room temperature before injection.
- Injections of **Stylage® XXL** are done with a needle/cannula according to the habits of the practitioner.
- Screw the needle/cannula firmly onto the Luer Lock hub of the syringe.
- Observe the aseptic and handling rules related to this type of intervention.
- Methods of injection:
 - . Slowly inject into the area indicated, using the sterile needle provided/the cannula chosen. It is recommended to use the technique of linear threading injection, multipoint injection, or a combination of both techniques. If an injection is too deep, the effectiveness of the correction will be reduced. If an injection is too superficial, discolouration of the skin or slight indurations may appear or the correction may be irregular.
- The amount to be injected depends on the skin defect to be corrected. Only the practitioner can define the dose to be injected for obtaining an optimal correction.
- Do not overcorrect. In the event of overcorrection, small indurations or irregular correction may appear.
- It is advised not to use more than 2.0 ml per treated area in a single session, and not to inject more than 20 ml of cross-linked hyaluronic acid per person per year.
- If the needle/cannula is blocked, do not increase the pressure required for injection. Replace the needle/cannula.
- After the injection, do not apply a cold pack; massage the treated area well in order to optimise the uniformity of the correction.

ADVICE TO PATIENTS

- Recommend using a sunscreen with a high protection factor during the two weeks following treatment.
- Advise the patient not to wear make-up during the 12 hours following the injection and avoid exposing the treated area to intense heat (UV, sauna, steam room) or extreme cold, at least until any potential post-injection swelling or redness disappear.

- The practitioner should tell the patient that he/she must keep him informed of any "abnormal" developments in the treated area (see adverse effects).
- Recommend that the patient avoid any massage or compression of the injected area for the first 3 days after the injection.
- Recommend that the patient avoid the sauna, steam room, and vigorous sports during the first week after injection.

ADVERSE EFFECTS

The practitioner must inform the patient of any potential immediate or delayed adverse effects that may occur following injection of the products, in particular (this list is not exhaustive):

- Inflammatory reactions such as redness, oedema, or erythema, potentially associated with itching and/or pain at the injection site, which usually resolve in less than a week.
- Haematomas.
- Induration or nodules, colouration or discolouration in the injected area.
- Poor efficacy or a weak filling effect.
- Local mobility of the implant.
- Rare cases of necrosis, abscesses, granuloma and hypersensitivity have been reported in literature following injections of hyaluronic acid. The patient must be informed about this.
- In patients who have a severe predisposition to allergies, dermatological disease, haemostasis disorder, or inflammatory disease, or in the event the precautions for use have not been observed, the incidence of adverse effects may increase.
- The patient must inform the practitioner of any adverse effects mentioned above that lasts for more than one week, or the appearance of any other adverse effect. The practitioner must report it to the reseller or manufacturer as soon as possible and should carry out an appropriate care.

STORAGE - EXPIRY DATE

- The products must be used before the expiration date indicated on the packaging and on the syringe.
- The products must be stored in their original packaging, at a temperature between 2°C and 25°C [36°F and 77°F], away from frost and light.
- If storage conditions are not observed, the product's performance can be compromised.

1. Remove the end cap of the syringe carefully.

2. Hold the ring at the end of the syringe and firmly insert the hub of the needle into the Luer-lock end of the syringe, by screwing until you feel some **strong counter pressure**.

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STJ1AM_B

OPDC341 Version 4 (11.2017)

This document is a photographic image of the product insert provided for your convenience by Med Supply Solutions. Med Supply Solutions does not claim ownership of the content contained within this product insert. The intellectual property rights, including but not limited to copyrights, trademarks, and any other applicable rights, pertaining to the content of this product insert, remain the sole property of the original manufacturer. Med Supply Solutions accepts no liability for any inaccuracies, omissions, or discrepancies that may be present in this document. Your use of this product insert is subject to our terms and conditions, which can be found on our website. By accessing and utilizing this document, you acknowledge and accept the terms and conditions outlined therein.

1-844-226-8277
info@medsupplysolutions.com
medsupplysolutions.com